



## Microalbuminuria (MAU) Diagnostic Kit



A low but abnormal level ( $> 30\text{mg/day}$ ) of MAU (urine micro-albumin) in urine is considered as an early marker of renal failure and cardiovascular complications in IDDM and NIDM. Persistent increase in urine albumin excretion rate between 30-300mg/day or 20-200mg/liter is the earliest clinically identifiable sign of risk of diabetic nephropathy and other vascular complications.

### Clinical Use of [MAU Test](#)

- ◆ Early marker for diabetic complications;
- ◆ Marker for renal failure in Hypertension;
- ◆ Related to cardiovascular diseases;
- ◆ The earliest signal and marker for renal failure.

### Parameters

Specimen	Urine
Specimen volume	80 $\mu\text{l}$
Cut off	20 $\mu\text{g/ml}$
Linear range	20-80 $\mu\text{g/ml}$
Precision	CV<15%
Accuracy	Relative deviation <20%
Interpretation time	5min
Test time	15s
Comparison	Compared with DiaSorin BN Prospec specific protein analyzer, linear correlation 0.92, agreement 94.2%

### Application

Endocrinology Dept, Nephrology Dept, Cardiology Dept, Clinical Laboratory, Health Examination Center

Tel: 86-21-68129666  
Fax: 86-21-68129222

Web site: [www.chemtronbio.com](http://www.chemtronbio.com)  
E-mail: [international@chemtronbio.com](mailto:international@chemtronbio.com)